



AZIENDA CERTIFICATA
SECONDO GLI STANDARD
DI QUALITÀ ISO 9001:2015

TECHNICAL DATASHEET

GLOVE TITS

Zootechny / Teat hygiene

Sealing gel for the dry phase

Features

High-power sealing filming gel designed to prevent environmental agents from penetrating into the teat canal of the quarter during dry and pre-birth period.

Application fields

GLOVE TITS is employed in the zootechny-veterinary sector.

How to use

Before use, thoroughly clean the nipples, disinfect the end and the papillary duct, then with paper or disposable wipes dry well. Using a disposable cup, immerse the last 3/4 cm of the nipple with the gel, wait 1 or 2 min. and then free the animal.

Concentration of use: as it is, without dilution

Allow to dry for: 1/2 min.

Features

Form: gel
density: 0,81 – 0,84 g/cm³ @25 °C

Packaging

Bottle: 0,5 kg
Bottle: 1 kg
Canister: 5 kg

Other information

GLOVE TITS meets the classification criteria of Regulation (EC) no. 1272/2008 (CLP). For specific information see safety data sheet.

GLOVE TITS is defined as dangerous for the different modes of transport according to the "United Nations Recommendations on the Transport of Dangerous Goods - Model Regulations". For specific information, see Section 14 of the Safety Data Sheet.

Attention

This technical data sheet is intended to be purely indicative. For the definition of the correct protocol for the use of Devidet s.r.l. always consult the Devidet s.r.l. Technical Office, expert technicians will be able to provide the most suitable solutions for specific needs. Product for professional use only. We decline all responsibility for damage resulting from improper or unauthorized use of the product. Before use, check that the product matches the type of application required. Before use, consult the safety data sheet or the information sheet.

Revisione 02, 23/02/2024

DEVIDET s.r.l.