

Protective gel for the dry and pre-birth period

*High-power protective
filming gel designed to
prevent environmental
agents from penetrating
the teat canal of the
quarter during dry and
pre-birth period.*

ADVANTAGES

- **UP TO 4/5 DAYS OF COVERAGE**
- **LOW COST PER COW**
- **WATERPROOF**
- **HIGH VISIBILITY**
- **1 LT. GLOVE TITS X100 TREATED COWS**

GUARANTEED BARRIER

APPLY BEFORE DRYING OFF THE COWS

Apply with a disposable cup. It dries very quickly forming a very resistant and highly durable elastomeric film, contributing to the reduction of the onset of mastitis.

HOW TO USE WITH THREE SIMPLE STEPS

- 1 - Check all the quarters, wearing disposable gloves, carefully clean the teats, including the bottom of the nipples.**
- 2 - Disinfect the bottom of the nipples and the papillary duct with paper or disposable wet wipes.**
- 3 - Using a disposable cup, immerse the last 3/4 cm of the nipple with **GLOVE TITS** wait 4/5 min. and then release the animal. For a better filming effect, a second immersion 20/30 seconds after the first treatment is recommended**



Conf. GR
500-1000-5000